

08/10/2020  
Thursday

## Expt I

### Determination of $\text{Na}_2\text{CO}_3$ and $\text{NaOH}$ in a mixture by titration.

#### Aim

To determine the amount of  $\text{Na}_2\text{CO}_3$  and  $\text{NaOH}$  in a mixture using  $\text{HCl}$  acid.

#### Principle

When a known volume of the mixture of  $\text{Na}_2\text{CO}_3$  and  $\text{NaOH}$  is treated with  $\text{HCl}$  using phenolphthalein indicator, at the end point all the hydroxide ions ( $\text{OH}^-$ ) and only half of the carbonate ions ( $\text{CO}_3^{2-}$ ) are reacted with  $\text{HCl}$  acid.

$A$  = all hydroxide ions + half of carbonate ions

$B$  = half the carbonate ions after phenolphthalein end point.

$2B$  = all carbonate ions

$A - B$  = all hydroxide ions

#### Procedure

##### Titration I: Standardization of $\text{HCl}$ .

Pipette out 20 ml of 0.1 N  $\text{Na}_2\text{CO}_3$  solution into a clean conical flask and add 2-3 drops of methyl orange indicator to the solution. Then titrate the solution against  $\text{HCl}$  acid taken in the burette. Record end point (burette reading) when colour changes from yellow to orange.

Repeat the titration till the concordant (2 consecutive burette readings exactly same) value is obtained.

Titration II : Estimation of  $\text{Na}_2\text{CO}_3$  and  $\text{NaOH}$  in a given mixture.

Dilute the given unknown solution to 100 ml in a std flask using distilled water. Pipette out 20 ml of this made up solution into a clean conical flask. Add 2-3 drops of phenolphthalein to the solution and titrate against standardized HCl. Record end point when disappearance of pink colour is observed. Let's consider this burette reading at 2 end points to be 'A' ml. To the same solution, add 2-3 drops of methyl orange and continue the titration till colour changes from yellow to orange, and note down burette reading as methyl orange end point. Consider it to be 'B' ml. Repeat the titration till the concordant value for A and B is obtained.

### Result

Amt of  $\text{Na}_2\text{CO}_3$  present in the given solution  
= 0.3089 g

Amt of  $\text{NaOH}$  present in the given solution  
= 0.0776 g

## TITRATION I :- Standardization of HCl.

Sl No.	Volume of $\text{Na}_2\text{CO}_3$ solution (ml)	Burette Reading		Concordant Values	Indica- tor
		Initial	Final		
1	20	0	19.5	19.5	Methyl orange.
2	20	0	19.5	19.5	

### Calculation

Volume of HCl = 19.5  $V_1$  ml (end point)  
Normality of HCl = 0.0512  $N_1$

Volume of  $\text{Na}_2\text{CO}_3 (V_2)$  = 20 mL  
Normality of  $\text{Na}_2\text{CO}_3 (N_2)$  = 0.05

$$\begin{aligned}\text{Normality of HCl } (N_1) &= \frac{N_2 V_2}{V_1} = \frac{20 \times 0.05}{19.5} \\ &= \cancel{0.0512} \quad \underline{\underline{0.0512 \text{ N}}}\end{aligned}$$

## TITRATION II - Estimation of $\text{Na}_2\text{CO}_3$ and $\text{NaOH}$ in a given mixture.

Sl No.	Volume of $\text{Na}_2\text{CO}_3$ unknown sol. (ml)	Burette Reading (ml)		
		Initial	Volume consumed for phenolphthalein end point (A)	Volume consumed for methyl orange end point after phenolphthalein end point (B)
1	20	0	19	30.4
2	20	0	19	30.4
CONCORDANT VALUE			(A) 19	(B) 11.4

### CALCULATION

#### I Estimation of amount of $\text{Na}_2\text{CO}_3$ .

$$\text{Volume of HCl } (V_1) = 2B = 2 \times 11.4 = 22.8 \text{ ml}$$

$$\text{Normality of HCl } (N_1) = N_1 = 0.0512$$

$$\text{Volume of mixture } (V_2) = 20 \text{ ml.}$$

$$\text{Normality of } \text{Na}_2\text{CO}_3 (N_2) = \frac{V_1 N_1}{V_2} = \frac{22.8 \times 0.0512}{20} = 0.0583 \text{ N}$$

$$\text{Amt of } \text{Na}_2\text{CO}_3 \text{ present in whole of given solution} = \frac{N_2 \times 53}{10} = \frac{0.0583 \times 53}{10} = 0.3089 \text{ g}$$

#### II Estimation of amt of $\text{NaOH}$

$$\text{Volume of HCl } (V_1) = (A - B) \text{ ml} = 7.6 \text{ ml.}$$

$$\text{Normality of HCl } (N_1) = 0.0512 \text{ N}$$

$$\text{Volume of mixture } (V_2) = 20 \text{ ml.}$$

$$\text{Normality of } \text{NaOH} (N_2) = \frac{V_1 \times N_1}{V_2} = \frac{7.6 \times 0.0512}{20} = 0.0194 \text{ N}$$

$$\text{Amt. of NaOH present in whole given solution} = \frac{N_2 \times 40}{10} = \frac{0.0194 \times 40}{10} = 0.0776 \text{ g}$$